



THE HEARING SOLUTION FOR YOUR INDIVIDUAL NEEDS.

Everyone's hearing is different. That is reflected in the extensive range of **in-the-ear hearing systems** Audio Service offers. Your hearing systems specialist can select the perfect solution for your personal hearing requirements from this comprehensive product range. You can experience the audio world of your chosen system for yourself during a non-binding trial period.

Simply consult your hearing system specialist. He or she will be happy to advise and support you.

You can make your own **personal hearing solution** a reality with Audio Service hearing systems. Many different factors demand consideration.

- Where exactly do your hearing deficits lie?
- Where and when do you want your hearing to be good?
- Which requirements arise from your life?

The **range of models** available and the **sophisticated technology** behind our hearing systems create the necessary preconditions for **good hearing**. Your hearing systems specialist is ideally qualified to advise you on the selection and customization process. Read about the technological performance that can be achieved with Audio Service in-the-ear hearing systems and the diverse range of possibilities they have to offer here.

Audio Service GmbH · Zeppelinstraße 9 · D-32051 Herford · Germany
www.audioservice.com

Art. No. 10536 809 / 05.2013 / Subject to technical modifications. Errors excepted. Colors may vary depending on the printing process.

DISCREET. INTELLIGENT. RELIABLE.

IN-THE-EAR HEARING SYSTEMS FROM AUDIO SERVICE

YOUR OWN PERSONAL HEARING SOLUTION.

HEARING SYSTEMS FROM AUDIO SERVICE

You want to enjoy good hearing again and to lead an active life – **in-the-ear hearing systems from Audio Service** provide reliable support in every situation. Those are excellent qualifications for secure, confident and long-lasting enjoyment of the world of hearing.

 **Audio Service**
hear · understand · communicate

www.audioservice.com



THE DISCREET TECHNICAL MARVEL IN THE AUDITORY CANAL

In-the-ear hearing systems provide a cosmetically discreet solution that will suit all seeking systems that are as inconspicuous as possible. In-the-ear hearing systems are worn directly in the auditory canal.

Hearing systems are produced to fit to your own specific auditory canal, which can be compared to a fingerprint. No two are alike. Great care is taken both in the production of your hearing system and its adjustment to match your personal hearing loss.

The technology used is all packed into the small in-the-ear hearing system. **A perfect fit, maximum discretion and custom configuration** enable a maximum of comfort for the wearer and the pleasure of unimpeded hearing.

A UNIQUE BLENDING OF SOUND THE CONCHA AND TECHNOLOGY

The acoustic quality familiar to every person is significantly influenced by the concha, which picks up sound coming from all directions, amplifies and concentrates it, and conveys it on to the eardrum. In-the-ear hearing systems retain this natural directional effect of the concha well, since they are located in the auditory canal.

Audio Service offers in-the-ear hearing systems in **a range of models and amplification versions** to suit different auditory canals, degrees of hearing loss and individual needs.

Six levels of technical sophistication ranging from entry-level through to innovative state-of-the-art technology provide solutions to match your personal requirements.

PLENTY OF TECHNOLOGY ON MINIMUM SPACE FOR INDIVIDUAL DEMANDS

In-the-ear hearing systems from Audio Service are **easy to use**. They **automatically** adapt to specific hearing situations and thus facilitate optimal comprehension of speech. The sophisticated technology used makes it possible to identify and automatically cancel out both background **noise** and any feedback that may occur. Program buttons allow you to **select from various hearing programs**.

A special **coating** ensures that the hearing systems are very easy to care for and keep clean. Special **filters** developed by Audio Service protect the delicate microphones and loudspeakers and keep the sensitive technology inside the hearing systems from soiling.

EXTRA SERVICE FOR YOUR EAR ACCESSORIES FOR MORE COMFORT

Selected in-the-ear hearing systems from Audio Service are optionally available with wireless functionality. Your hearing systems can then be networked with other audio devices via a **Bluetooth remote control** system. Audio signals transmitted directly to your hearing systems – cable-free and in stereo quality – open up a new dimension of state-of-the-art infotainment.

Optionally, a **remote control unit** can be supplied with various in-the-ear hearing systems for easy adjustment of hearing programs and volume.

Discover
the world of
hearing on the
Internet:
audioservice.com